

AllerFlow | Gluten

Hygiena™ AllerFlow Gluten is a rapid and convenient test for detecting gluten residue on surfaces as part of an allergen control program. Combining Hygiena's convenient sample collection swab design with classic lateral flow technology, AllerFlow Gluten makes gluten testing easier than ever.


Benefits:

- Convenient and quick sample collection requires minimal training
- Pre-moistened swab enables reliable collection to ensure accurate results
- Simple two-part kit includes all required materials - no pipettes, bottles, scales, or heating blocks required
- Pre-measured, self-contained extraction buffer provides consistent results with minimal handling
- Easy-to-interpret lateral flow test results in 10 minutes
- Specific to the toxic fraction of gluten from wheat, barley, rye, and oat
- No cross-reactivity with soy, rice, or corn
- Up to 15-month shelf life at room or refrigerated temperature (2 – 25 °C)

Sensitivity:

- Sensitive to 5 µg of gluten per swab within 10 minutes

Positive Result:



Negative Result:



Catalog No.

Description

Quantity

ALF-GL


AllerFlow Gluten Kit

25



Find support documents, instructional videos, and more at www.hygiena.com

AllerFlow Gluten Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

Related Products



EnSURE™ Monitoring System

EnSURE™ is a high-sensitivity luminometer that measures multiple quality tests such as ATP, TVC, Coliform, *E. coli*, and more, all on one system. The EnSURE Monitoring System includes SureTrend Data Analysis Software.



AllerSnap™ High-Sensitivity Allergen Prevention Test

AllerSnap™ is a colorimetric allergen prevention test for broad range allergen screening in 15-30 minutes, making it ideal for everyday use.



SuperSnap™ High-Sensitivity Surface ATP Test

SuperSnap™ is Hygiena's most sensitive ATP surface test. Designed to work with all Hygiena luminometers, SuperSnap is an all-in-one test which detects extremely low levels of ATP, enabling its use in food and beverage processing facilities as an allergen cross-contamination prevention tool.

Catalog No.	Description	Quantity
ALS-100	AllerSnap™ High-Sensitivity Allergen Prevention Test	100
ENSURE	EnSURE™ Monitoring System with SureTrend Software	1
SUS3000/SUS3000X	SuperSnap™ High-Sensitivity Surface ATP Test	100



Find support documents, instructional videos, and more at www.hygiena.com