

Features & Benefits

- Advanced photodiode technology and unique chemistry for maximum sensitivity and durability
- Intuitive smart navigation makes it easy for anyone to use
- Functions like a smartphone
- Syncs data to SureTrend Cloud
- Supports multiple languages



- Responsive 5-inch shatter-proof touch screen
- Train remote teams with built-in screen sharing capability
- Built-in tilt sensor for more accurate readings
- Auto shut-off for smart battery management
- Simulate competitor upper and lower RLU limits

Compatible Products

UltraSnap™ | Surface ATP Test



Catalog No. US2020 | Qty: 100 [Instructional Video Online](#)

SuperSnap™ | High-Sensitivity Surface ATP Test



Catalog No. SUS3000 | Qty: 100
Catalog No. SUS3000X | Qty: 100 [Instructional Video Online](#)

AquaSnap TOTAL™ | Water ATP Test



Catalog No. AQ-100X | Qty: 100 [Instructional Video Online](#)

AquaSnap FREE™ | Water ATP Test



Catalog No. AQ-100FX | Qty: 100 [Instructional Video Online](#)

CalCheck™ | LED Calibration Verification Device



Catalog No. CAL | Qty: 1 [Instructional Video Online](#)

MicroSnap™ | Indicator Organism Tests



Catalog No. MS2-EB | Qty: 100
Catalog No. MS2-ECOLI | Qty: 100
Catalog No. MS2-TOTAL | Qty: 100
Catalog No. MS2-COLIFORM | Qty: 100 [Instructional Video Online](#)

ZymoSnap ALP™ | ALP Pasteurization Verification Test



Catalog No. ZS-ALP-100 | Qty: 100 [Instructional Video Online](#)

CROSS CHECK™ | ACP Cooking Efficiency Test



Catalog No. CX-3000 | Qty: 100

Product Name	Catalog No.	Quantity
EnSURE™ Touch	ETOUCH	1
UltraSnap™ Surface ATP Test	US2020	100
SuperSnap™ High-Sensitivity Surface ATP Test	SUS3000	100
SuperSnap™ High-Sensitivity Surface ATP Test	SUS3000X	100
AquaSnap™ Total Water ATP Test	AQ-100X	100
AquaSnap™ Free Water ATP Test	AQ-100FX	100
MicroSnap™ Coliform & E. coli Enrichment Device	MS1-CEC	100
MicroSnap™ Enhanced Nutrient Broth	MS1-N-BROTH-9ML	100
MicroSnap™ Coliform Detection Device	MS2-COLIFORM	100
MicroSnap™ E. coli Detection Device	MS2-ECOLI	100
MicroSnap™ Enterobacteriaceae Enrichment Device	MS1-EB	100

Product Name	Catalog No.	Quantity
MicroSnap™ Enterobacteriaceae Detection Device	MS2-EB	100
MicroSnap™ Total Enrichment Device	MS1-TOTAL	100
MicroSnap™ Total Detection Device	MS2-TOTAL	100
ZymoSnap ALP Alkaline Phosphatase Test	ZS-ALP-100	100
CrossCheck Acid Phosphatase Test	CX-3000	100
CalCheck LED Calibration Verification Device	CAL	1
ATP Positive Control Kit	CK25	25
ZymoSnap ALP Positive Control Kit	ZS-ALP-PC	5
Small Format Incubator	INCUBATOR	1
Lab Format Incubator	INCUBATOR2	1
BAX System Q7	ASY2018	1
BAX System X5	MCH2007	1

EnSURE Touch is factory calibrated and runs a calibration check on each startup. Recommended calibration control kits available.

help.hygiene.com

[YouTube.com/HygieneTV](https://www.youtube.com/HygieneTV)

[Vimeo.com/HygieneTV](https://www.vimeo.com/HygieneTV)

Biogenetics DIAGNOSTICS S.R.L. • Viale Germania 1/B • 35020 Ponte San Nicolò PD • Tel. 049 8968186
diagnostics@biogenetics.it • www.biogenetics.it

Americas:

Hygiene Headquarters

941 Avenida Acaso
Camarillo, CA 93012
Phone: 1-805-388-8007
Fax: 1-805-388-5531
info@hygiene.com

Hygiene Canada

2650 Meadowvale Blvd Unit 14
Mississauga, Ontario L5N 6M5
Phone: 1-833-494-4362 (Toll-free)
or 1-416-686-7962

Hygiene Mexico, S.A. de C.V.

Calle 3 Anegas 409 Bodega 5, Col. Nueva
Industrial Vallejo, Delegación Gustavo A.
Madero, C.P. 07700, CDMX, México.
Phone: +52-5281-4108 or +52-5281-4146

International:

Hygiene International

8, Woodshots Meadow
Watford, Hertfordshire
WD18 8YU, UK
Phone: +44 (0)1923-818821
Fax: +44 (0)1923-818825

Hygiene (Shanghai) Trading Co., Ltd.

Rm.7K, No.518, Shangcheng Rd.
Pudong New District
Shanghai, China 200120
Phone: +86 21-5132-1081, +86 21-5132-1077,
+86 21-5132-1078
Fax: +86 21-5132-1080

Hygiene Diagnóstica España S.L.

P. I. Parque Plata, Calle Cañada Real 31-35,
41900, Camas, Sevilla, Spain
Phone: +34 954-08-1276

EnSURE™ Touch

Sanitation Verification.
Designed for You.



hygiene.com/ensuretough-fb





EnSURE™ Touch

Today's food safety professionals face significant challenges: complex supply chains, increased audit and inspection scrutiny, and the ever-present risk of contamination and recall. Professionals need to ensure processes are free of contamination and efficiently perform product safety interventions.

EnSURE™ Touch is the next-generation quality monitoring system. It provides rapid and accurate sanitation verification to support your FSMA program. Featuring a 5-inch touch screen, wireless sync technology, and cloud-based software, the EnSURE Touch is designed to adapt to your workplace and provides the data you need for complex multi-location operations, audit and risk management, and recall prevention.

EnSURE Touch Can Help You

- Verify sanitation of surfaces, liquids, and product testing in as little as 10 seconds
- Test and compare the effectiveness of sanitation products and procedures
- Train and monitor sanitation crew and QC/QA employees
- Support SSOP and objectively prove cleaning standards have been met
- Quickly identify when replaceable parts like belts and brushes should be replaced
- Identify problem areas in your facility that need more attention
- Provide on-the-spot training and feedback

In Order To:

- Meet FSMA, GFSI, and ISO standards
- Protect your brand reputation
- Avoid financially damaging lawsuits
- Reduce risk of food poisoning, facility closure, and product recalls
- Lower facility testing costs by bringing testing in-house
- Provide evidence of due diligence - HACCP



Innovation at Your Fingertips

Quick test: Run a test in as few as two touches

Locations: Create testing points with customized fields and pass/fail limits on the go without needing desktop software

Plans: Easily group testing locations into plans and randomize test locations to ensure sampling coverage

Users: Manage user access to differentiate personnel from managers. Enhanced administrator control

Results & Reports: Easily view testing results and trend analysis

Search: Quickly find specific testing plans and locations on the fly, limiting time between tests

Sync: Wirelessly synchronize data to SureTrend Cloud

Calibration: Factory-calibrated and preset checks to stay GMP/ISO-compliant


Retest: Swipe through failed test locations and instantly run re-tests to show effectiveness of corrective action

MicroSnap™: Test and show results for organisms like Coliform, *E. coli*, Enterobacteriaceae, and total viable count*



Enzymes: Test and show results for enzymes like Alkaline Phosphatase and Acid Phosphatase*

Settings: Customize settings for language, custom fields, test types, and simulate competitor RLU limits

Specifications: Size: 7" x 3" x 1" • Weight: 0.7 lbs • Memory: 2GB internal + cloud storage • Connectivity: USB or Wi-Fi*
Language: Supports multiple languages • Battery: Rechargeable Lithium-Ion battery with USB-C charging *



* Battery life varies based on use and configuration • *MicroSnap and Enzyme functionality coming soon in future software update • *USB connectivity coming soon in future software update



The data you need to master risk management all in one place.

SureTrend Cloud is an easy-to-use data analysis software that comes free with every EnSURE Touch. It enables users to monitor testing results and manage user access across multiple facilities, making risk management easier than ever.

Features & Benefits:

- Communicate department performance and achievement to your team and to management
- Fully interactive dashboard allows for quick data viewing, showing you everything at a glance
- Easily manage one or hundreds of instruments from a single SureTrend Cloud account
- Monitor testing across one or multiple facilities
- Full online support and resource center only one click away
- Track re-tests to discover opportunities for retraining or identify equipment that may need replacing
- Quickly create custom reports to trend testing performance
- Automatic updates to SureTrend Cloud ensure you're always running the latest version

suretrend.hygiena.com

